

Operation Manual

Shodex OHpak SB-807 HQ

(Please read this manual carefully before using the column to keep its life and performance well.)

1. Introduction

Shodex OHpak SB-807 HQ is a gel filtration chromatography (GFC) column for water-soluble huge polymers. It is possible to measure a water-soluble huge polymer without partial exclusion even if its molecular weight distribution of polymer-side reaches several ten million or more.

2. Instructions in handling <Important>

Caution!

* Take notice of keeping instructions about the solvents and the reagents used with the column not to occur problems related to losing your health or leaking.

Attention!

* Use the column within the regular range of flow rate, pressure and temperature. There is a danger of deteriorating the column suddenly when it is handled beyond the permissible range even for a short time. See the clause "Usable conditions" about the permissible range.

3. Specifications

Column size	: 8.0 mm I.D. × 300 mm L.
Column material (inside)	: SUS 316
Packing material	: Poly(hydroxymethacrylate)-type porous particle
In-column solvent (initial)	: Ion exchanged water
Exclusion limit M.W.	: ca. 500,000,000 (Estimated value from the extrapolated calibration curve by standard pullulans.)
Number of theoretical plates	: $\geq 1,500$ per column

4. Usable conditions

Flow rate	: ≤ 1.5 ml/min * 0.5~1.0ml/min is suitable usually.
Pressure	: ≤ 0.5 MPa per column
Temperature	: 4~60 °C
pH	: 3~10 * The eluent should be within the range of pH6~10 to avoid the corrosion of the device and the column when chloride ion is included.
Eluent	: Aqueous solution of salt or buffer solution is usually used. <Typical salt> Sodium chloride, sodium nitrate, sodium sulfate, potassium sulfate, ammonium sulfate. <Typical buffer> Phosphate buffer, Tris-HCl buffer, acetate buffer, citrate buffer. * Acetonitrile or methanol can be added within 30 vol% to suppress the adsorption of the hydrophobic sample.

Attention!

- 1) Do not remove the end fittings of the column under any circumstances.
- 2) Do not make a strong impact on the column: such as hitting or dropping on the floor.
- 3) Replace the solvent in the chromatograph with the eluent to be used before connecting the column.
- 4) Connect the column so that the flow direction corresponds to the arrow mark on the tag.
- 5) When the column is not used for a month or more, replace the in-column solvent with ion exchanged water, close each end with a stopper, and store it at room temperature.
- 6) Borate buffer is not desirable because there is a possibility of forming the complex with the hydrophilic surface of the packing material.
- 7) Confirm there is no cloudiness in the solution when the aqueous solution of salt is used added with acetonitrile or methanol.
- 8) Adjust the total concentration of the salts to 0.5M or less. In general, the range of 0.05~0.3M is suitable.
- 9) Install a guard column (SB-807G) upstream to protect the main column from deterioration.